

April 23, 2020

### Additional Information Regarding Cleaning Agents

S.L.P.'s products should be wiped with a non-corrosive to plastic cleanser before use.

The sensors are not to be autoclaved or soaked in disinfectants.

S.L.P. has tested and can recommend on the following agents for the cleaning of its sensors:

Product Name	Brand Name	Active Ingredients
SaniZide Plus <sup>®</sup> Germicidal Wipes	Safetec	<ul style="list-style-type: none"> <li>• Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.105%</li> <li>• Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.105%</li> </ul>
CaviWipes <sup>™</sup> Disinfecting Towelettes	Metrex	<ul style="list-style-type: none"> <li>• Diisobutylphenoxyethoxyethyl dimethyl benzyl ammonium chloride 0.28%</li> <li>• Isopropanol 17.20%</li> </ul>
Super Sani-Cloth <sup>®</sup> Germicidal Disposable Wipe	PDI	<ul style="list-style-type: none"> <li>• n-alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.25%</li> <li>• n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.25%</li> </ul>
Protex <sup>®</sup> Ultra Disinfectant Wipes	Parker	<ul style="list-style-type: none"> <li>• Octyl decyl dimethyl ammonium chloride 0.0909%</li> <li>• Dioctyl dimethyl ammonium chloride 0.0364%</li> <li>• Didecyl dimethyl ammonium chloride 0.0545%</li> <li>• Alkyl (C14, 50%; C12, 40%; C16, 10%) dimethyl benzyl ammonium chloride 0.1212%</li> </ul>

Opti-Cide <sup>3</sup> ® Surface Wipes	Bristol	<ul style="list-style-type: none"> <li>• n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride 0.154%</li> <li>• n-Alkyl (68% C12, 32% C14) dimethyl ethyl benzyl ammonium chloride 0.154%</li> <li>• Isopropanol 21.000%</li> </ul>
Sani-Cloth® AF3 Germicidal Disposable Wipe	PDI	<ul style="list-style-type: none"> <li>• n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chlorides 0.14%</li> <li>• n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chlorides 0.14%</li> </ul>
Meliseptol® Wipes sensitive	B. Braun	<ul style="list-style-type: none"> <li>• Propan-1-ol 17%</li> <li>• Didecyldimethylammoniumchloride 0.23%</li> </ul>
Rubbing Alcohol		<ul style="list-style-type: none"> <li>• Isopropyl Alcohol 70%</li> <li>• Ethyl Alcohol (Ethanol) 70%</li> </ul>

Additionally, the following agents were tested on S.L.P.'s belts. These agents to be used according to the manufacturer instructions. After using these agents, the belts can be machine washed with household laundry detergent (e.g. Tide®, Seventh Generation®, etc.).

Product Name	Brand Name	Active Ingredients
Professional Surface Disinfectant Spray	Purell	Ethyl Alcohol 29.4%
Disinfectant Spray Plus	Mediclean	<ul style="list-style-type: none"> <li>• o-Phenylphenol 0.22%</li> <li>• Diisobutylphenoxyethoxy ethyl dimethyl benzyl ammonium chloride monohydrate 0.70%</li> </ul>



Inbal Raphael  
QA/RA Manager  
S.L.P. Ltd.