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
**Sensores de Ronquidos**

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**Schnarchsensoren**

Thank you for purchasing a SleepSense® sensor.

### Description

Snoring sensors respond to snoring and other sounds in the audio range picked up through the skin, and convert them to a small analog voltage that provides a clear, reliable indication of the presence of these sounds. Piezo or dynamic sensors do not require a battery or interface. Electret microphones are powered by the system. They are easily activated by plugging the output connector into the appropriate input on the head box.

- Dynamic snore sensors offer a wide frequency range, and sensitivity to other respiratory sounds besides snoring.
- Piezo sensors are robust, and relatively insensitive to ambient noise.
- Electret microphones offer a wide frequency response and high sensitivity.

### Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

### Applying the Sensor

- The sensor is attached to the patient's neck.
- To position the sensor, place the sensor with the opening side or the protruding side (for the Piezo Sensor), against the skin and use medical grade tape (intended for facial skin) to hold the sensor in place. Make sure the sensor is tight against the skin to ensure that vibrations are clearly transmitted to the sensor during snoring.
- Ask the patient to simulate snoring, and check that clear, strong signals are being transmitted before leaving the patient and starting the study.

### Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.
- Do not soak the sensor in disinfectants.

### Technical Specifications:

Sensor:	Dynamic / Piezo	Electret
Cable:	White insulation wire	
Sensor design:	Round flexible TPE disc	
Recommended filter settings:	High pass: 50 HZ / Low pass: 500 HZ or as recommended in the system's instructions	
Signal Output:	Approximately 1 mV for loud snoring	Approximately 50 mV for loud snoring
Operating conditions:	5°C (40°F) - 40°C (104°F)	
Storage temperature:	-20°C (-4°F) - 60°C (140°F)	
Operating and storage humidity:	5% - 95% (Non-condensing)	

### Warnings and Precautions

- **Caution:** US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense sensors are for professional use only.
- SleepSense sensors may be used only in conjunction with an approved recording system.
- SleepSense sensors are only intended for use on healthy skin.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.

#### Limited Warranty

SLP warrants the snoring sensor to be free of defects in materials and workmanship for a period of six months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our Distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This Warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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