Inductive Sensors

Italiano

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Specifications are subject to change without notice

S.L.P. Inc.


Inductive Sensors

SleepSense®

The Sensor Solution Company

User Manual
Multi-Rip (Two Channels + Sum)
Manuel de l'utilisateur
Multi-Rip (Deux Canaux + Somme)
Manual del Usuario
Multi-Rip (Dos Canales + Suma)
Manuale per l'Utente
Multi-Rip (Due Canali + Somma)
Bedienungsanleitung
Multi-Rip (Zwei Kanäle + Summe)

Thank you for purchasing a SleepSense® sensor.

Description

The inductance of the inductive bands is proportional to the cross section area of the band's closed loop around the body. It therefore changes as the patient breathes. The electronics in the interface box measure the thoracic and abdominal effort signals and their sum (sum model) to provide accurate and reliable respiratory effort signals used to assist the technician in identifying desynchronized chest and abdominal breathing.

Instructions for Use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

Applying the Sensor

The patient should lay on his/her back for optimal positioning of both bands.
For use with SLP inductive bands.

Positioning the Sensor (See Figure A)

1. Strap the appropriately sized sensor bands around the patient's chest and abdomen, over any clothing.
2. Attach the first wire to the snap that is on the same side of the band as the Velcro® patch.
3. While holding this end to the left of the patient's centerline, wrap the bands around the patient's chest and abdomen. For smaller patients move the starting point further to the left.
4. Pull the other end over the first end and secure to the band on the right side of the centerline. The band should be snug but not too tight.
5. Attach the second wire to the exposed snap. Secure both wires about 3" away from the connection.

- **When removing the sensor band do not pull on the wires. Remove from the snap.**

Interface Box Setup (See Figure B)

- **The chest interface cable MUST be connected to chest inputs to activate the unit. Once connected, a green light appears indicating the device is ON.**
 - A red light appears when the interface box is operating and the battery power is running low. When this occurs, it is recommended to replace the batteries within the next few days.
1. Plug the black connector of the chest interface cable into the designated socket on the interface box. Align the notch of the plug with the grove of the socket. When you hear a click, the plug is properly connected.
 2. Plug the black connector of the abdomen interface cable into the designated socket on the interface box. Align the notch of the plug with the grove of the socket. When you hear a click, the plug is properly connected.
 3. Plug the black output connector into the designated socket on the interface box. Align the notch of the plug with the grove of the socket. When you hear a click, the plug is properly connected.
 - Connect the output pins to the corresponding inputs in the system. The black pin is the common reference for all channels.
- **When removing the interface cables do not pull on the wires. Remove from the connector.**
 - **Check that clear, strong signals are being transmitted before leaving the patient.**

Cleaning the Sensor

- The bands are machine washable.
- Wipe the Interface Cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the bands and Interface Cable are thoroughly dry before reusing it.
- Do not soak the Interface Cable in disinfectants.

Technical Specifications:

Items #:	Box 9201, Chest cable PS9203, Abdomen Cable PS9202, Output cable PS9204
Recommended filter settings:	High pass: 0.1 Hz / Low pass: 40 Hz or as recommended in the systems' instructions
Operating conditions:	5°C (40°F) - 40°C (104°F)
Storage temperature:	-20°C (-4°F) - 60°C (140°F)
Operating and storage humidity:	5% - 95% (Non-condensing)

! Warnings and Precautions
• Caution: US Federal Law restricts this device to sale, distribution and use by or on the order of a physician. • SleepSense sensors are for professional use only. • SleepSense sensors may be used only in conjunction with an approved recording system. • SleepSense sensors are only intended for use on healthy skin. • Discontinue use if sensor shows signs of wear, damage, or exposed metal. Attention: contains a lithium battery - dispose of properly.

Limited Warranty

SLP warrants the Multi-Rip box and cables to be free of defects in materials and workmanship for a period of twelve months from the date purchased. The sole liability of SLP and our Distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our Distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This Warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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