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
Thank you for purchasing a SleepSense® sensor.

Description

Comprising of a single heat sensitive element, the oral respiration-flow sensor is attached onto the pressure cannula under the patient's nose. The sensor produces a clear, reliable signal that is directly proportional to the temperature changes of respiratory oral air flow. This sensor may be connected to the respiration flow channel on most recorders.

Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system. They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.

Applying the Sensor

The Respiration Flow Sensor is attached onto the pressure cannula which is placed on the patient's upper lip. The sensors are delicate and must be handled with care.

Positioning the Sensor

1. Bend a 1/4" section of the end of the black-cupped wire upwards by 90 degrees.
Fold them forward and down over the cannula body.
Lock the sensor to the cannula by pressing on the wires.
2. Align the tip of the oral prong with the center of the mouth/upper lip.
3. Draw the wires on either side of the sensor over the patient's ears and use the clasp under the chin to secure the wires in place.

Recommendations:

- The sensor prong should not come into contact with the skin. Contact with the skin will significantly attenuate the amplitude of the signal.
- Tape the wire on both cheeks to reduce strain and enhance reliability.
- Check that clear, strong signals are being transmitted before leaving the patient.

Cleaning the Sensor

- Wipe the sensor and cable with a non-corrosive (to plastic) cleanser to clean before use.
- Make sure the complete sensor assembly is thoroughly dry before reusing it.
- To sterilize the sensors, use standard gas sterilization procedures.
- Do not autoclave or soak the sensor in disinfectants.

Technical Specifications:

For a complete list of technical specifications please contact us or visit our website.

Cable:	long white PVC insulation lead wire
Recommended filter settings:	High pass: 0.1 Hz / Low pass: 70 Hz or as recommended in the system's instructions
Operating conditions:	5°C (40°F) - 40°C (104°F)
Storage temperature:	-20°C (-4°F) - 60°C (140°F)
Operating and storage humidity:	5% - 95% (Non-condensing)

! Warnings and Precautions
• SleepSense sensors are for professional use only. • Caution: US Federal Law restricts this device to sale, distribution and use by or on the order of a physician. • SleepSense sensors may be used only in conjunction with an approved recording system. • SleepSense sensors are only intended for use on healthy skin. • Discontinue use if sensor shows signs of wear, damage, or exposed metal.

Limited Warranty
SLP warrants the respiration-flow sensor to be free of defects in materials and workmanship for a period of six months from the date purchased. The sole liability of SLP and our distributor(s) is limited to replacement or repair of the product at the option of SLP with no charge for parts or labor if any part is proven to be defective in workmanship, performance, or materials during the warranty period. Under no circumstances shall SLP or our distributor(s) be liable for any loss of revenues or damage, direct, consequential, or incidental, including loss of profit, property damage, or personal injury arising from the use of, or the inability to use this product. This warranty is intended only for the original buyer and is in lieu of all other warranties or previous agreements, expressed or implied. This warranty is rendered void if the product is used for other than its intended purpose or is subject to abuse, misuse, tampering, neglect, or unauthorized modifications. Use of this product constitutes acceptance of this warranty in total.

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