

# SleepSense®

The Sensor Solutions Company

ThermoCan™ 7-Ft Child Thermal Cannula  
ThermoCan™ Canule Thermique, Enfant, 210cm  
ThermoCan™ Cánula Térmica, Niño, 210cm  
ThermoCan™ Cannula Termica, Bambino, 210cm  
ThermoCan™ Thermal Kanüle, Kind, 210cm



Use by: **LOT** 00000000  
YYYY / MM **REF** 15814

**S.L.P. Ltd.**

62 Mordechai Anilevitz St.  
Tel-Aviv, 6706016, Israel  
Tel: +972 3 6870640  
Fax: +972 3 5371282

**S.L.P. Inc.**

2511 Technology Dr.,  
Suite 112 Elgin, IL 60124, USA  
Tel: +1 (630) 513 7479  
Fax: +1 (630) 513 8478  
E-mail: info@sleepsense.com

**EC REP**

**Arazy Group GmbH**  
The Squaire 12,  
Am Flughafen,  
60549 Frankfurt am Main,  
Germany

**CH REP**

**Arazy Group  
Switzerland GmbH**  
Bruderholzallee 53, 4059  
Basel, Switzerland

**UK**

UK Responsible Person:


**Medes Ltd.**  
5 Beaumont Gate,  
Shenley Hill, Herts,  
Radlett D77AR, England,  
United Kingdom

Made in China

Item N15814 - ver.3.1 [side A]

[www.sleepsense.com](http://www.sleepsense.com)

Toll-free(U.S.): 888-slp-sens (757 7367)



### Directions for Use

- Connect the thermal cannula using SleepSense child interface cable by inserting metal plugs on either side of the nose piece. (*figure 1*)
- Position cannula on face by inserting cannula tips into nostrils. (*figure 2*)
- Connect one tube to a Pressure Sensor and the other to a Capnograph.
- Connect interface cable to the PSG system.
- Ask the patient to breath through his nose and mouth and ensure that clear signals are recorded.

### Indications for use

SleepSense sleep-lab sensors provide a qualitative measure of sounds, respiratory-effort, flow, body position or limb movements, for recording on an approved data acquisition system.

They are intended for use on children and adult patients who are screened during sleep disorder studies at a sleep laboratory or the patient's home.



### Warnings and Precautions

- **Caution:** This product is for diagnostic purposes only and is not to be used as an apnea monitor or in a life sustaining situation.  
US Federal Law restricts this device to sale, distribution and use by or on the order of a physician.
- SleepSense sensors are for professional use only.
- SleepSense sensors may be used only in conjunction with an approved recording system.
- Discontinue use if sensor shows signs of wear, damage, or exposed metal.
- SleepSense sensors are only intended for use on healthy skin.